

BUSA® High Speed Burs






Stryker® 3.16mm TPS™ Shank Replacements



Compatible with High Speed Drill Systems: Stryker® CORE & TPS™

Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

STRYKER® | Compatible with High Speed Drill Systems: Stryker® CORE & TPS™
CONMED® LINVATEC® | Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II







Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
KM-M8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Stryker® 5120-107-30, 5220-107-30, 5420-107-30, ConMed® Linvatec® 7021-237, Medtronic® Midas Rex® M-8				
				
KM-AM8U	Match Head	mm	2.9	3.8
		in	0.11	0.15
Replaces Stryker® 5130-107-30, 5230-107-30, 5430-107-30, ConMed® Linvatec® 7021-337, Medtronic® Midas Rex® AM-8				
				
KM-M51U	Metal Cutter (Diamond Wheel)	mm	25.4	0.8
		in	1.00	0.03
Replaces Stryker® 5120-130-250, ConMed® Linvatec® 7021-292, Medtronic® Midas Rex® M-51				
				
KM-M52U	Metal Cutter (Carbide)	mm	3.1	17.5
		in	0.12	0.69
Replaces Stryker® 5120-80-30, 5420-80-30, ConMed® Linvatec® 7021-294, Medtronic® Midas Rex® M-52				
				
KM-U295	Oval	mm	4.0	8.0
		in	0.16	0.31
Replaces Stryker® 5120-15-140, ConMed® Linvatec® 7021-295, Medtronic® Midas Rex® N/A				
				

BUSA® High Speed Burs

Stryker® 3.16mm TPS™ Shank Replacements

Compatible with High Speed Drill Systems: Stryker® CORE & TPS™

Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
KM5110-15-40	Oval	mm	4.0	9.5
		in	0.16	0.37
Replaces Stryker® 5110-15-40, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				
KM-M340DU	Round Diamond	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5120-12-40, ConMed® Linvatec® 7021-264, Medtronic® Midas Rex® M-340D				
				
KM-M330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5120-10-30, 5220-10-30, ConMed® Linvatec® 7021-226, Medtronic® Midas Rex® M-330				
				
KM5120-10-230	Round Fluted Extended Length	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5120-10-230, 5220-10-230, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				
KM-M340U	Round Fluted	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5120-10-40, 5220-10-40, ConMed® Linvatec® 7021-227, Medtronic® Midas Rex® M-340				
				
KM5120-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5120-10-240, 5220-10-240, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				

STRYKER® | Compatible with High Speed Drill Systems: Stryker® CORE & TPS™

CONMED® LINVATEC® | Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

BUSA® High Speed Burs







Stryker® 3.16mm TPS™ Shank Replacements



Compatible with High Speed Drill Systems: Stryker® CORE & TPS™

Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

STRYKER® | Compatible with High Speed Drill Systems: Stryker® CORE & TPS™
CONMED® LINVATEC® | Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II


Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
KM-M35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5120-10-50, 5220-10-50, 5420-10-50, ConMed® Linvatec® 7021-228, Medtronic® Midas Rex® M-35				
				
KM-M33U	Round Fluted	mm	6.0	6.0
		in	0.24	0.24
Replaces Stryker® 5120-10-60, 5220-10-60, ConMed® Linvatec® 7021-229, Medtronic® Midas Rex® M-33				
				
KM-AM330U	Round Fluted	mm	3.0	3.0
		in	0.12	0.12
Replaces Stryker® 5130-10-30, 5230-10-30, ConMed® Linvatec® 7021-326, Medtronic® Midas Rex® AM-330				
				
KM5130-10-240	Round Fluted Extended Length	mm	4.0	4.0
		in	0.16	0.16
Replaces Stryker® 5130-10-240, 5230-10-240, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				
KM-AM35U	Round Fluted	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5130-10-50, 5230-10-50, 5430-10-50, ConMed® Linvatec® 7021-328, Medtronic® Midas Rex® AM-35				
				
KM5130-10-250	Round Fluted Extended Length	mm	5.0	5.0
		in	0.20	0.20
Replaces Stryker® 5130-10-250, 5230-10-250, ConMed® Linvatec® N/A, Medtronic® Midas Rex® N/A				
				

BUSA® High Speed Burs

Stryker® 3.16mm TPS™ Shank Replacements

Compatible with High Speed Drill Systems: Stryker® CORE & TPS™

Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

Stryker®				
Order Number	Description	Unit mm/inch	Head Diameter	Head Length
KM-AM32U	Round Fluted	mm	7.6	7.6
		in	0.30	0.30
Replaces Stryker® 5130-10-75, 5230-10-75, ConMed® Linvatec® 7021-331, Medtronic® Midas Rex® AM-32				
				







STRYKER® | Compatible with High Speed Drill Systems: Stryker® CORE & TPS™

CONMED® LINVATEC® | Compatible with ConMed® Linvatec® UltraPower®; Medtronic® Midas Rex® Classic I & II

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-C0901	Oval - Carbide	mm	44.5	4.0	7.9	8
		in	1.75	0.16	0.31	
Replaces Stryker® 277-10-62, 277-010-062						
						1 per pack, Sterile
14-C0905	Oval - Carbide	mm	44.5	5.5	10.3	10
		in	1.75	0.22	0.41	
Replaces Stryker® N/A						
						1 per pack, Sterile
16-C0901	Oval - Carbide	mm	54.5	4.0	7.9	8
		in	2.15	0.16	0.31	
Replaces Stryker® N/A						
						1 per pack, Sterile
14-C0104	Round - Carbide	mm	44.5	1.8	1.8	6
		in	1.75	0.07	0.07	
Replaces Stryker® 277-10-418, 277-010-418						
						1 per pack, Sterile
14-C0105	Round - Carbide	mm	44.5	2.0	2.0	6
		in	1.75	0.08	0.08	
Replaces Stryker® 277-10-225, 277-010-225						
						1 per pack, Sterile
14-C0106	Round - Carbide	mm	44.5	2.3	2.3	8
		in	1.75	0.09	0.09	
Replaces Stryker® 277-10-423, 277-010-423						
						8 Shannon 1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-C0107	Round - Carbide	mm	44.5	3.0	3.0	6
		in	1.75	0.12	0.12	
Replaces Stryker® 277-10-230, 277-010-230						
						1 per pack, Sterile
14-C0108	Round - Carbide	mm	44.5	4.0	4.0	8
		in	1.75	0.16	0.16	
Replaces Stryker® 277-10-240, 277-010-240						
						1 per pack, Sterile
14-C0109	Round - Carbide	mm	44.5	5.0	5.0	8
		in	1.75	0.20	0.20	
Replaces Stryker® 277-10-250, 277-010-250						
						1 per pack, Sterile
14-C0110	Round - Carbide	mm	44.5	6.0	6.0	8
		in	1.75	0.24	0.24	
Replaces Stryker® N/A						
						1 per pack, Sterile
16-C0107	Round - Carbide	mm	54.5	3.0	3.0	6
		in	2.15	0.12	0.12	
Replaces Stryker® N/A						
						1 per pack, Sterile
16-C0108	Round - Carbide	mm	54.5	4.0	4.0	8
		in	2.15	0.16	0.16	
Replaces Stryker® N/A						
						1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-C0308	Straight Cross Cut - Carbide	mm	44.5	1.6	4.5	6
		in	1.75	0.06	0.17	
Replaces Stryker® 1607-2-103, 1607-002-103						
						560 1 per pack, Sterile
14-C0405	Straight Cross Cut - Carbide	mm	44.5	1.6	5.9	6
		in	1.75	0.06	0.23	
Replaces Stryker® N/A						
						1 per pack, Sterile
14-C0606	Straight Cross Cut - Carbide	mm	44.5	1.2	4.1	6
		in	1.75	0.05	0.16	
Replaces Stryker® 277-10-212, 277-010-212						
						1 per pack, Sterile
14-C0608	Tapered Cross Cut - Carbide	mm	44.5	1.6	4.3	6
		in	1.75	0.06	0.17	
Replaces Stryker® 277-10-216, 277-010-216						
						1 per pack, Sterile
14-C0616	Tapered Cross Cut - Carbide	mm	44.5	2.1	4.8	6
		in	1.75	0.08	0.19	
Replaces Stryker® 277-10-221, 277-010-221						
						1 per pack, Sterile
16-D0130	Round - Diamond	mm	54.5	2.0	2.0	Fine
		in	2.15	0.08	0.08	
Replaces Stryker® 1608-6-87, 1608-006-087						
						1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-D0131	Round - Diamond	mm	54.5	2.4	2.4	Fine
		in	2.15	0.09	0.09	
Replaces Stryker® 1608-6-89, 1608-006-089						
						1 per pack, Sterile
16-D0132	Round - Diamond	mm	54.5	3.0	3.0	Medium
		in	2.15	0.11	0.11	
Replaces Stryker® 1608-6-91, 1608-006-091						
						1 per pack, Sterile
16-D0134	Round - Diamond	mm	54.5	3.5	3.5	Medium
		in	2.15	0.14	0.14	
Replaces Stryker® 1608-6-93, 1608-006-093						
						1 per pack, Sterile
16-D0135	Round - Diamond	mm	54.5	4.0	4.0	Medium
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-6-95, 1608-006-095						
						1 per pack, Sterile
16-D0137	Round - Diamond	mm	54.5	5.0	5.0	Medium
		in	2.15	0.20	0.20	
Replaces Stryker® 1608-6-99, 1608-006-099						
						1 per pack, Sterile
14-S0402	Cylinder Stainless Steel	mm	44.5	6.1	11.1	14
		in	1.75	0.24	0.44	
Replaces Stryker® 1607-2-1, 1607-002-001						
						1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S0402	Cylinder Stainless Steel	mm	54.5	6.1	11.1	14
		in	2.15	0.24	0.44	
Replaces Stryker® 1608-2-1, 1608-002-001						
						1 per pack, Sterile
15-S1101	Flame - Stainless Steel	mm	48.0	6.0	12.7	14
		in	1.89	0.24	0.50	
Replaces Stryker® 1607-2-47, 1607-002-047						
						1 per pack, Sterile
14-S0201	Inverted Cone Stainless Steel	mm	44.5	3.1	3.3	10
		in	1.75	0.12	0.13	
Replaces Stryker® 1607-2-57, 1607-002-057						
						1 per pack, Sterile
16-S1805	Lindemann Stainless Steel	mm	54.5	1.6	9.0	4
		in	2.15	0.06	0.35	
Replaces Stryker® N/A						
						1 per pack, Sterile
16-S1801	Lindemann Stainless Steel	mm	54.5	2.3	9.5	4
		in	2.15	0.09	0.37	
Replaces Stryker® 1608-2-21, 1608-002-021						
						1 per pack, Sterile
19-S1802	Lindemann Stainless Steel	mm	66.7	2.2	24.3	4
		in	2.63	0.09	0.96	
Replaces Stryker® 1608-2-43, 1608-002-043						
						1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
23-S1809	Lindemann Stainless Steel	mm	70.0	2.3	22.0	4
		in	2.76	0.09	0.87	
Replaces Stryker® N/A						
						1 per pack, Sterile
22-S1808	Lindemann Stainless Steel	mm	75.0	2.3	35.0	4
		in	2.95	0.09	1.38	
Replaces Stryker® N/A						
						1 per pack, Sterile
14-S0903	Oval - Stainless Steel	mm	44.5	4.0	9.5	10
		in	1.75	0.16	0.37	
Replaces Stryker® 1607-2-35, 1607-002-035						
						1 per pack, Sterile
14-S0909	Oval - Stainless Steel	mm	44.5	6.0	9.5	14
		in	1.75	0.24	0.37	
Replaces Stryker® 1607-2-3, 1607-002-003						
						1 per pack, Sterile
16-S0903	Oval - Stainless Steel	mm	54.5	4.0	9.5	10
		in	2.15	0.16	0.37	
Replaces Stryker® 1608-2-35, 1608-002-035						
						1 per pack, Sterile
16-S0909	Oval - Stainless Steel	mm	54.5	6.0	9.5	14
		in	2.15	0.24	0.37	
Replaces Stryker® 1608-2-3, 1608-002-003						
						1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
14-S1201	Pear - Stainless Steel	mm	44.5	6.1	10.0	14
		in	1.75	0.24	0.39	
Replaces Stryker® 1607-2-5, 1607-002-005						
						1 per pack, Sterile
16-S1201	Pear - Stainless Steel	mm	54.5	6.1	10.0	14
		in	2.15	0.24	0.39	
Replaces Stryker® 1608-2-5, 1608-002-005						
						1 per pack, Sterile
16-S0106	Round - Stainless Steel	mm	54.5	2.3	2.3	6
		in	2.15	0.09	0.09	
Replaces Stryker® 1608-6-147, 1608-006-147						
						1 per pack, Sterile
16-S0107	Round - Stainless Steel	mm	54.5	3.2	3.2	6
		in	2.15	0.126	0.126	
Replaces Stryker® 1608-6-149, 1608-006-149						
						1 per pack, Sterile
16-S0108	Round - Stainless Steel	mm	54.5	4.0	4.0	8
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-6-137, 1608-006-137						
						1 per pack, Sterile
16-S0109	Round - Stainless Steel	mm	54.5	5.1	5.1	8
		in	2.15	0.20	0.20	
Replaces Stryker® 1608-6-141, 1608-006-141						
						1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S0110	Round - Stainless Steel	mm	54.5	6.1	6.1	8
		in	2.15	0.24	0.24	
Replaces Stryker® 1608-6-143, 1608-006-143						
						1 per pack, Sterile
16-S0115	Round - Stainless Steel	mm	54.5	2.4	2.4	8
		in	2.15	0.095	0.095	
Replaces Stryker® 1608-2-9, 1608-002-009						
						1 per pack, Sterile
16-S0116	Round - Stainless Steel	mm	54.5	3.0	3.0	20
		in	2.15	0.12	0.12	
Replaces Stryker® 1608-6-163, 1608-006-163						
						1 per pack, Sterile
16-S0117	Round - Stainless Steel	mm	54.5	3.2	3.2	8
		in	2.15	0.126	0.126	
Replaces Stryker® 1608-2-11, 1608-002-011						
						1 per pack, Sterile
16-S0118	Round - Stainless Steel	mm	54.5	4.0	4.0	24
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-6-165, 1608-006-165						
						1 per pack, Sterile
16-S0119	Round - Stainless Steel	mm	54.5	4.0	4.0	6
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-6-151, 1608-006-151						
						1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S0120	Round - Stainless Steel	mm	54.5	4.0	4.0	12
		in	2.15	0.16	0.16	
Replaces Stryker® 1608-2-55, 1608-002-055						
 <p style="text-align: right;">1 per pack, Sterile</p>						
16-S0121	Round - Stainless Steel	mm	54.5	4.6	4.6	6
		in	2.15	0.18	0.18	
Replaces Stryker® 1608-6-153, 1608-006-153						
 <p style="text-align: right;">1 per pack, Sterile</p>						
16-S0122	Round - Stainless Steel	mm	54.5	4.6	4.6	8
		in	2.15	0.18	0.18	
Replaces Stryker® 1608-6-139, 1608-006-139						
 <p style="text-align: right;">1 per pack, Sterile</p>						
16-S0123	Round - Stainless Steel	mm	54.5	4.5	4.5	12
		in	2.15	0.177	0.177	
Replaces Stryker® 1608-2-13, 1608-002-013						
 <p style="text-align: right;">1 per pack, Sterile</p>						
16-S0126	Round - Stainless Steel	mm	54.5	7.1	7.1	8
		in	2.15	0.28	0.28	
Replaces Stryker® 1608-6-145, 1608-006-145						
 <p style="text-align: right;">1 per pack, Sterile</p>						
16-S1401	Shannon Style Tapered - Stainless Steel	mm	54.5	2.1	11.1	3
		in	2.15	0.08	0.44	
Replaces Stryker® 1608-2-101, 1608-002-101						
 <p style="text-align: right;">1 per pack, Sterile</p>						

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Standard & Oral Surgery

Stryker® J-Notch Replacements

Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
16-S1601	Swanson Reamer Stainless Steel	mm	56.5	2.0	10.5	5
		in	2.23	0.08	0.41	
Replaces Stryker® 1608-2-131, 1608-002-131						
						1 per pack, Sterile
17-S1602	Swanson Reamer Stainless Steel	mm	62.0	3.0	13.5	5
		in	2.44	0.12	0.53	
Replaces Stryker® 1608-2-133, 1608-002-133						
						1 per pack, Sterile
17-S1603	Swanson Reamer Stainless Steel	mm	62.0	4.0	13.5	6
		in	2.44	0.16	0.53	
Replaces Stryker® 1608-2-135, 1608-002-135						
						1 per pack, Sterile
18-S1701	Wire Pass Drill Stainless Steel	mm	64.5	1.0	10.0	2
		in	2.54	0.039	0.39	
Replaces Stryker® 1608-2-63, 1608-002-063						
						1 per pack, Sterile
18-S1703	Wire Pass Drill Stainless Steel	mm	65.5	1.5	6.9	2
		in	2.57	0.059	0.27	
Replaces Stryker® 1608-2-59, 1608-002-059						
						1 per pack, Sterile
19-S1705	Wire Pass Drill Stainless Steel	mm	66.7	2.0	6.4	2
		in	2.63	0.079	0.25	
Replaces Stryker® 1608-2-61, 1608-002-061						
						1 per pack, Sterile

STRYKER® | Compatible with Stryker® Micro Drill Systems: RemB™, CORE, TPS™, F1, TPX & Command 2

BUSA® Micro Burs - Oral Surgery

Stryker® Flat Shank Replacements







Compatible with Stryker® Impaction Drill Systems - 100K, CMD, CORE & Oral Max

Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KMS958N	Round - Carbide	mm	58.4	2.3	2.3	8
		in	2.30	0.09	0.09	
Replaces Stryker® 2296-101-423						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)
KMS94701N	Tapered Fissure Cross Cut - Carbide	mm	44.5	1.2	4.2	6
		in	1.75	0.05	0.17	
Replaces Stryker® 2296-100-212						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)
KMS94702N	Tapered Fissure Cross Cut - Carbide	mm	44.5	1.6	4.4	6
		in	1.75	0.06	0.17	
Replaces Stryker® 2296-100-216						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)
KMS94703N	Tapered Fissure Cross Cut - Carbide	mm	44.5	2.1	4.8	6
		in	1.75	0.08	0.19	
Replaces Stryker® 2296-100-221						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)
KMS95701N	Tapered Fissure Cross Cut - Carbide	mm	58.4	1.2	4.2	6
		in	2.30	0.05	0.17	
Replaces Stryker® 2296-101-212						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)
KMS95702N	Tapered Fissure Cross Cut - Carbide	mm	58.4	1.6	4.4	6
		in	2.30	0.06	0.17	
Replaces Stryker® N/A						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)

STRYKER® | Compatible with Stryker® Impaction Drill Systems - 100K, CMD, CORE & Oral Max

BUSA® Micro Burs - Oral Surgery Stryker® Flat Shank Replacements

Compatible with Stryker® Impaction Drill Systems - 100K, CMD, CORE & Oral Max







Stryker®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
KMS94703LN	Tapered Fissure Long Cross Cut - Carbide	mm	44.5	2.1	7.5	6
		in	1.75	0.08	0.30	
Replaces Stryker® 2296-100-321						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)
KMS95702LN	Tapered Fissure Long Cross Cut - Carbide	mm	58.4	1.6	6.0	6
		in	2.30	0.06	0.24	
Replaces Stryker® 2296-101-316						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)
KMS941703LN	Tapered Fluted Carbide	mm	44.5	2.1	7.0	8
		in	1.75	0.08	0.27	
Replaces Stryker® 2296-100-421						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)
KMS951703LN	Tapered Fluted Carbide	mm	58.4	2.1	7.0	8
		in	2.30	0.08	0.27	
Replaces Stryker® 2296-101-421						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)
KM94901N	Oval - Carbide	mm	48.0	4.0	8.0	8
		in	1.89	0.16	0.31	
Replaces Stryker® 2296-100-26, 2296-100-026						
						Non-Sterile
KMS95901N	Oval - Carbide	mm	58.4	4.0	8.0	8
		in	2.30	0.16	0.31	
Replaces Stryker® 2296-101-26, 2296-101-026						
						5 per pack, Sterile (Individual Sterile Packs with Lot Numbers)

STRYKER® | Compatible with Stryker® Impaction Drill Systems - 100K, CMD, CORE & Oral Max

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
03-C0901	Oval - Carbide	mm	48.0	4.0	8.0	8
		in	1.89	0.16	0.31	
Replaces ConMed® Linvatec® 5091-36, 5091-136, 5091-236, M40B						
						
03-C0905	Oval - Carbide	mm	48.0	5.5	10.2	10
		in	1.89	0.22	0.40	
Replaces ConMed® Linvatec® 5091-38, 5091-138, 5091-238, M55B						
						
09-C0901	Oval - Carbide	mm	70.0	4.0	8.0	8
		in	2.76	0.16	0.31	
Replaces ConMed® Linvatec® 5092-36, 5092-136, 5092-236, L40B						
						
09-C0905	Oval - Carbide	mm	70.0	5.5	10.2	10
		in	2.76	0.22	0.40	
Replaces ConMed® Linvatec® 5092-38, 5092-138, 5092-238, L55B						
						
12-C0901	Oval - Carbide	mm	95.0	4.0	8.0	8
		in	3.74	0.16	0.31	
Replaces ConMed® Linvatec® 5093-36, 5093-136, 5093-236, X40B						
						
12-C0905	Oval - Carbide	mm	95.0	5.5	10.2	10
		in	3.74	0.22	0.40	
Replaces ConMed® Linvatec® 5093-38, 5093-138, 5093-238, X55B						
						







CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems







ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0101	Round - Carbide	mm	51.0	1.0	1.0	6
		in	2.01	0.04	0.04	
Replaces ConMed® Linvatec® 5091-20, 5091-120, 5091-220						
						
04-C0103	Round - Carbide	mm	51.0	1.5	1.5	6
		in	2.01	0.06	0.06	
Replaces ConMed® Linvatec® 5091-22, 5091-122, 5091-222, M15R						
						
04-C0104	Round - Carbide	mm	51.0	1.8	1.8	6
		in	2.01	0.07	0.07	
Replaces ConMed® Linvatec® 5091-21, 5091-121, 5091-221, M18R						
						
04-C0105	Round - Carbide	mm	51.0	2.0	2.0	6
		in	2.01	0.08	0.08	
Replaces ConMed® Linvatec® 5091-24, 5091-124, 5091-224, M20R						
						
04-C0106	Round - Carbide	mm	51.0	2.3	2.3	8
		in	2.01	0.09	0.09	
Replaces ConMed® Linvatec® 5091-23, 5091-123, 5091-223, M23R, 5091-25, 5091-125, 5091-225, M24R						
						
04-C0107	Round - Carbide	mm	51.0	3.0	3.0	8
		in	2.01	0.12	0.12	
Replaces ConMed® Linvatec® 5091-26, 5091-126, 5091-226, M30R						
						

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

CONMED® LINVATEC®

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems







ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0108	Round - Carbide	mm	51.0	4.0	4.0	8
		in	2.01	0.16	0.16	
Replaces ConMed® Linvatec® 5091-28, 5091-128, 5091-228, M40R						
						
04-C0109	Round - Carbide	mm	51.0	5.0	5.0	10
		in	2.01	0.20	0.20	
Replaces ConMed® Linvatec® 5091-30, 5091-130, 5091-230, M50R						
						
04-C0110	Round - Carbide	mm	51.0	6.0	6.0	8
		in	2.01	0.24	0.24	
Replaces ConMed® Linvatec® M60R						
						
04-C0170	Round - Carbide	mm	51.0	6.5	6.5	10
		in	2.01	0.26	0.26	
Replaces ConMed® Linvatec® 5091-32, 5091-132, 5091-232, M65R						
						
09-C0109	Round - Carbide	mm	70.0	5.0	5.0	10
		in	2.76	0.20	0.20	
Replaces ConMed® Linvatec® 5092-30, 5092-130, 5092-230, L50R						
						
09-C0110	Round - Carbide	mm	70.0	6.0	6.0	8
		in	2.76	0.24	0.24	
Replaces ConMed® Linvatec® L60R						
						

CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems







ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-C0170	Round - Carbide	mm	70.0	6.5	6.5	10
		in	2.76	0.26	0.26	
Replaces ConMed® Linvatec® 5092-32, 5092-132, 5092-232, L65R						
						
10-C0101	Round - Carbide	mm	71.5	1.0	1.0	6
		in	2.81	0.04	0.04	
Replaces ConMed® Linvatec® 5092-20, 5092-120, 5092-220, L10R						
						
10-C0103	Round - Carbide	mm	72.0	1.5	1.5	6
		in	2.83	0.06	0.06	
Replaces ConMed® Linvatec® 5092-22, 5092-122, 5092-222, L15R						
						
10-C0105	Round - Carbide	mm	72.5	2.0	2.0	6
		in	2.85	0.08	0.08	
Replaces ConMed® Linvatec® 5092-24, 5092-124, 5092-224, L20R						
						
10-C0107	Round - Carbide	mm	73.0	3.0	3.0	8
		in	2.87	0.12	0.12	
Replaces ConMed® Linvatec® 5092-26, 5092-126, 5092-226, L30R						
						
10-C0108	Round - Carbide	mm	73.0	4.0	4.0	8
		in	2.87	0.16	0.16	
Replaces ConMed® Linvatec® 5092-28, 5092-128, 5092-228, L40R						
						

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

CONMED® LINVATEC®

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems







ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
12-C0105	Round - Carbide	mm	95.0	2.0	2.0	6
		in	3.74	0.08	0.08	
Replaces ConMed® Linvatec® 5093-24, 5093-124, 5093-224, X20R						
						
12-C0108	Round - Carbide	mm	95.0	4.0	4.0	8
		in	3.74	0.16	0.16	
Replaces ConMed® Linvatec® 5093-28, 5093-128, 5093-228, X40R						
						
12-C0109	Round - Carbide	mm	95.0	5.0	5.0	10
		in	3.74	0.20	0.20	
Replaces ConMed® Linvatec® 5093-30, 5093-130, 5093-230, X50R						
						
04-C0601	Tapered - Carbide	mm	51.0	0.9	5.2	6
		in	2.01	0.035	0.20	
Replaces ConMed® Linvatec® 5091-15, 5091-115, 5091-215						
						
04-C0606	Tapered - Carbide	mm	51.0	1.2	4.2	6
		in	2.01	0.05	0.17	
Replaces ConMed® Linvatec® M12T						
						
04-C0623	Tapered - Carbide	mm	51.0	1.2	6.0	6
		in	2.01	0.05	0.24	
Replaces ConMed® Linvatec® N/A						
						

CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-C0607	Tapered - Carbide	mm	51.0	1.6	4.3	6
		in	2.01	0.06	0.17	
Replaces ConMed® Linvatec® 5091-18, 5091-118, 5091-218						
						
04-C0609	Tapered - Carbide	mm	51.0	1.5	10.0	6
		in	2.01	0.06	0.39	
Replaces ConMed® Linvatec® 5091-02, 5091-102, 5091-202, M15T						
						
04-C0612	Tapered Cross Cut - Carbide	mm	51.0	2.1	4.9	6
		in	2.01	0.08	0.19	
Replaces ConMed® Linvatec® 5091-13, 5091-113, 5091-213, M18T						
						
04-C0615	Tapered - Carbide	mm	51.0	2.1	4.8	6
		in	2.01	0.08	0.19	
Replaces ConMed® Linvatec® N/A						
						
04-C0617	Tapered - Carbide	mm	51.0	2.0	7.5	6
		in	2.01	0.08	0.30	
Replaces ConMed® Linvatec® 5091-06, 5091-106, 5091-206, M27T						
						
06-C0619	Tapered - Carbide	mm	57.0	2.0	19.0	6
		in	2.24	0.08	0.75	
Replaces ConMed® Linvatec® 5091-14, 5091-114, 5091-214, M20T						
						







Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

CONMED® LINVATEC®

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
10-C0617	Tapered - Carbide	mm	73.0	2.0	7.5	6
		in	2.87	0.08	0.30	
Replaces ConMed® Linvatec® 5092-06, 5092-106, 5092-206, L27T						
						
11-C0619	Tapered - Carbide	mm	80.0	2.0	19.0	6
		in	3.15	0.08	0.75	
Replaces ConMed® Linvatec® 5092-14, 5092-114, 5092-214, L20T						
						
03-C0312	Straight Cross Cut - Carbide	mm	48.0	4.0	10.0	6
		in	1.89	0.16	0.39	
Replaces ConMed® Linvatec® 5091-54, 5091-154, 5091-254, M40S						
						
04-C0302	Straight Cross Cut - Carbide	mm	51.0	1.0	5.0	6
		in	2.01	0.04	0.20	
Replaces ConMed® Linvatec® 5091-01, 5091-101, 5091-201, M10S						
						
04-C0306	Straight Cross Cut - Carbide	mm	51.0	1.5	7.5	6
		in	2.01	0.06	0.30	
Replaces ConMed® Linvatec® 5091-04, 5091-104, 5091-204, M15S						
						
04-C0307	Straight - Carbide	mm	51.0	1.6	4.3	6
		in	2.01	0.06	0.17	
Replaces ConMed® Linvatec® 5091-17, 5091-117, 5091-217, M17S						
						







CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery

2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems







ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
10-C0306	Straight - Carbide	mm	73.0	1.5	7.5	6
		in	2.87	0.06	0.30	
Replaces ConMed® Linvatec® 5092-04, 5092-104, 5092-204, L13S, L17S						
						
10-C0310	Straight Cross Cut - Carbide	mm	73.0	2.1	10.0	8
		in	2.87	0.08	0.39	
Replaces ConMed® Linvatec® 5092-52, 5092-152, 5092-252, L20S						
						
03-D0135	Round - Diamond	mm	48.0	4.0	4.0	Medium
		in	1.89	0.16	0.16	
Replaces ConMed® Linvatec® 5091-264, 5091-164, 5091-264, M40D						
						
04-D0139	Round - Diamond	mm	51.0	2.0	2.0	Medium
		in	2.01	0.08	0.08	
Replaces ConMed® Linvatec® 5091-60, 5091-160, 5091-260, M20D						
						
04-D0133	Round - Diamond	mm	51.0	3.0	3.0	Medium
		in	2.01	0.12	0.12	
Replaces ConMed® Linvatec® 5091-62, 5091-162, 5091-262						
						
09-D0135	Round - Diamond	mm	70.0	4.0	4.0	Medium
		in	2.76	0.16	0.16	
Replaces ConMed® Linvatec® 5092-64, 5092-164, 5092-264, L40D, L60D						
						

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

CONMED® LINVATEC®

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems







ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-D0141	Round - Diamond	mm	70.0	5.0	5.0	Coarse
		in	2.76	0.20	0.20	
Replaces ConMed® Linvatec® 5092-66, 5092-166, 5092-266, L50D						
						
09-D0138	Round - Diamond	mm	70.0	6.0	6.0	Coarse
		in	2.76	0.24	0.24	
Replaces ConMed® Linvatec® 5092-68, 5092-168, 5092-268						
						
10-D0139	Round - Diamond	mm	73.0	2.0	2.0	Medium
		in	2.87	0.08	0.08	
Replaces ConMed® Linvatec® 5092-60, 5092-160, 5092-260, L20D						
						
10-D0133	Round - Diamond	mm	73.0	3.0	3.0	Medium
		in	2.87	0.12	0.12	
Replaces ConMed® Linvatec® 5092-62, 5092-162, 5092-262, L30D						
						
08-S1802	Lindemann Stainless Steel	mm	67.0	2.3	25.0	4
		in	2.64	0.09	0.98	
Replaces ConMed® Linvatec® 5091-12, 5091-112, 5091-212, L17L						
						
09-S0115	Round - Stainless Steel	mm	70.0	2.4	2.4	8
		in	2.76	0.091	0.091	
Replaces ConMed® Linvatec® N/A						
						

CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems




ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
09-S0120	Round - Stainless Steel	mm	70.0	4.0	4.0	12
		in	2.76	0.16	0.16	
Replaces ConMed® Linvatec® N/A						
						
12-S0153	Round - Stainless Steel	mm	95.0	3.1	3.1	10
		in	3.74	0.12	0.12	
Replaces ConMed® Linvatec® N/A						
						
04-S1501	Shannon Style - Tapered Stainless Steel	mm	51.0	2.3	12.0	4
		in	2.01	0.09	0.47	
Replaces ConMed® Linvatec® N/A						
						
06-S1601	Swanson Reamer Stainless Steel	mm	57.0	2.0	10.0	5
		in	2.24	0.08	0.39	
Replaces ConMed® Linvatec® 5091-42, 5091-142, 5091-242, M20F						
						
07-S1602	Swanson Reamer Stainless Steel	mm	62.0	3.0	13.5	5
		in	2.44	0.12	0.53	
Replaces ConMed® Linvatec® 5091-43, 5091-143, 5091-243, M30F						
						
07-S1605	Swanson Reamer Stainless Steel	mm	62.0	4.0	16.5	6
		in	2.44	0.16	0.65	
Replaces ConMed® Linvatec® 5091-44, 5091-144, 5091-244, M40F						
						

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

CONMED® LINVATEC®

BUSA® Small Bone/Micro Burs Standard & Oral Surgery 2.35 mm Straight Shank - ConMed® Linvatec® Replacements

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems

ConMed® Linvatec®						
Order Number	Description	Unit	Total Length	Head Diameter	Head Length	No. of Flutes Grit Size
04-S1701	Wire Pass Drill Stainless Steel	mm	51.0	1.0	19.0	2
		in	2.01	0.04	0.75	
Replaces ConMed® Linvatec® 5091-47, 5091-147, 5091-247, M10W						
						
04-S1703	Wire Pass Drill Stainless Steel	mm	51.0	1.5	19.0	2
		in	2.01	0.059	0.75	
Replaces ConMed® Linvatec® 5091-48, 5091-148, 5091-248, M15W						
						
04-S1704	Wire Pass Drill Stainless Steel	mm	51.0	2.0	19.0	2
		in	2.01	0.079	0.75	
Replaces ConMed® Linvatec® 5091-49, 5091-149, 5091-249, M20W						
						

CONMED® LINVATEC®

Compatible with Small Bone/Micro Drill Systems: ConMed® Linvatec® - MicroPower®, MicroChoice®, Micro 100™, MicroE® & Surgairtome Two®; BUSA® & MicroAire® Small Bone/Micro Drill Systems