



NEW-MED s.r.l.

Piazza Santa Maria, 18
26855 Lodi Vecchio (LO) - Italy

+39 0371752374






info@new-med.eu
<https://www.new-med.eu/>

P.iva/C.F: IT-12502690964

Rea: LO-2666383

Large Bone Saw Blades - Oscillating/Sagittal BR2000™ (Aggressive Cut) MEDICAL - Proprietary (K5 Series - DePuy® Synthes® and Zimmer Biomet®)

Compatible with DePuy® Synthes® Oscillating Saw Systems: Air Power Line, Battery Power Line II, & Compact Air Drive II and Trauma Recon; Zimmer Biomet® Oscillating Saw System: Universal

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
K5-1290-35	mm	90.0	12.5	0.89	0.89
	in	3.54	0.49	0.035	0.035
Substitute for DePuy® Synthes® 05.002.201S, 05.003.100S, Zimmer Biomet® 13-0900-89A-F1					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					
K5-1290-39	mm	90.0	12.5	0.99	0.99
	in	3.54	0.49	0.039	0.039
Substitute for DePuy® Synthes® 05.002.405S, 05.003.102S, Zimmer Biomet® 13-0901-00A-F1					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					
K5-1990-35	mm	90.0	19	0.89	0.89
	in	3.54	0.75	0.035	0.035
Substitute for DePuy® Synthes® 05.002.001S, 05.003.110S, Zimmer Biomet® 19-0900-89A-F1					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					
K5-2590-35	mm	90.0	25	0.89	0.89
	in	3.54	0.98	0.035	0.035
Substitute for DePuy® Synthes® 05.002.101S, 05.003.130S, Zimmer Biomet® 25-0900-89A-F1					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					
K5-2590-58	mm	90.0	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
Substitute for DePuy® Synthes® 05.002.105S, 05.003.139S, Zimmer Biomet® 25-0901-47A-B1					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					