



NEW-MED s.r.l.

Piazza Santa Maria, 18  
26855 Lodi Vecchio (LO) - Italy

+39 0371752374



info@new-med.eu  
<https://www.new-med.eu/>

P.iva/C.F: IT-12502690964

Rea: LO-2666383


## Large Bone Saw Blade - Reciprocating Replacement - Zimmer Biomet®

Compatible with Zimmer Biomet® Reciprocating Saw System: Universal

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
ZB12-0761-00U-R1	mm	70.0	12.5	0.61	0.90
	in	2.76	0.49	0.024	0.035
<i>Replaces 12-0761-00U-R1</i>					
					<i>Tooth Offset</i>  <i>Sterile</i> <i>Double Sided</i>
ZB12-0901-19U-R1	mm	89.4	12.2	0.90	1.14
	in	3.52	0.48	0.035	0.045
<i>Replaces 12-0901-19U-R1</i>					
					<i>Tooth Offset</i>  <i>Sterile</i> <i>Single Sided</i>






## Sternum Saw Blade - Proprietary

Compatible with ConMed® Linvatec® Sagittal Saw Systems: MicroPower™, MicroChoice®, Micro 100™ & Micro E®, MicroAire® - Small Bone Sagittal Saw System; Brasseler USA® - All Small Bone Sagittal Saw System

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
KM-3231SS	mm	41.9	19.0	0.40	0.71
	in	1.65	0.75	0.016	0.028
<i>Replaces N/A</i>					
					<i>Sterile</i>

# Large Bone Saw Blades - Oscillating/Sagittal BR2000™ (Aggressive Cut) - Proprietary (K5 Series - DePuy® Synthes® and Zimmer Biomet®)

Compatible with DePuy® Synthes® Oscillating Saw Systems: Air Power Line, Battery Power Line II, & Compact Air Drive II and Trauma Recon; Zimmer Biomet® Oscillating Saw System: Universal

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>K5-1290-35</b>	mm	90.0	12.5	0.89	0.89
	in	3.54	0.49	0.035	0.035
<i>Substitute for DePuy® Synthes® 05.002.201S, 05.003.100S, Zimmer Biomet® 13-0900-89A-F1</i>					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					
<b>K5-1290-39</b>	mm	90.0	12.5	0.99	0.99
	in	3.54	0.49	0.039	0.039
<i>Substitute for DePuy® Synthes® 05.002.405S, 05.003.102S, Zimmer Biomet® 13-0901-00A-F1</i>					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					
<b>K5-1990-35</b>	mm	90.0	19	0.89	0.89
	in	3.54	0.75	0.035	0.035
<i>Substitute for DePuy® Synthes® 05.002.001S, 05.003.110S, Zimmer Biomet® 19-0900-89A-F1</i>					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					
<b>K5-2590-35</b>	mm	90.0	25	0.89	0.89
	in	3.54	0.98	0.035	0.035
<i>Substitute for DePuy® Synthes® 05.002.101S, 05.003.130S, Zimmer Biomet® 25-0900-89A-F1</i>					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					
<b>K5-2590-58</b>	mm	90.0	25	1.47	1.47
	in	3.54	0.98	0.058	0.058
<i>Substitute for DePuy® Synthes® 05.002.105S, 05.003.139S, Zimmer Biomet® 25-0901-47A-B1</i>					
 <p style="text-align: right;">No Tooth Offset</p> <p style="text-align: right;">Sterile</p>					