









# BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements

**Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2**

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM3-51</b>	mm	18.5	6.5	0.38	0.43
	in	0.73	0.26	0.015	0.017
Replaces Stryker® 2296-3-51, 2296-003-051					
					Thin
<b>KM3-101</b>	mm	10	9	0.38	0.43
	in	0.39	0.35	0.015	0.017
Replaces Stryker® 2296-3-101, 2296-003-101					
					Thin
<b>KM3-102</b>	mm	18.5	16.5	0.38	0.43
	in	0.73	0.65	0.015	0.017
Replaces Stryker® 2296-3-102, 2296-003-102					
					Thin
<b>KM3-103</b>	mm	15	7	0.38	0.43
	in	0.59	0.28	0.015	0.017
Replaces Stryker® 2296-3-103, 2296-003-103					
					Thin
<b>KM3-104</b>	mm	34.5	16.5	0.38	0.43
	in	1.36	0.65	0.015	0.017
Replaces Stryker® 2296-3-104, 2296-003-104					
					Thin
<b>KM3-105</b>	mm	18.5	9	0.38	0.43
	in	0.73	0.35	0.015	0.017
Replaces Stryker® 2296-3-105, 2296-003-105					
					Thin

Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2

**STRYKER®**

# BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements









**STRYKER®**

Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2

**Compatible with Stryker® Oscillating & Sagittal Saw Systems:**






**SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2**

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM3-106</b>	mm	34.5	13	0.38	0.43
	in	1.36	0.51	0.015	0.017
Replaces Stryker® 2296-3-106, 2296-003-106					
					Thin
<b>KM3-108</b>	mm	39	13	0.38	0.43
	in	1.54	0.51	0.015	0.017
Replaces Stryker® 2296-3-108, 2296-003-108					
					Thin
<b>KM3-111</b>	mm	25	9	0.38	0.43
	in	0.98	0.35	0.015	0.017
Replaces Stryker® 2296-3-111, 2296-003-111					
					Thin
<b>KM3-114</b>	mm	18.5	7	0.38	0.43
	in	0.73	0.28	0.015	0.017
Replaces Stryker® 2296-3-114, 2296-003-114					
					Thin
<b>KM3-115</b>	mm	29.5	7	0.38	0.43
	in	1.16	0.28	0.015	0.017
Replaces Stryker® 2296-3-115, 2296-003-115					
					Thin
<b>KM3-125</b>	mm	31	9	0.38	0.43
	in	1.22	0.35	0.015	0.017
Replaces Stryker® 2296-3-125, 2296-003-125					
					Thin



# BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements

Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
 SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM3-206</b>	mm	34.5	12	0.38	0.43
	in	1.36	0.47	0.015	0.017
Replaces Stryker® 2296-3-206, 2296-003-206					
					Thin
<b>KM3-212</b>	mm	18	5.5	0.38	0.43
	in	0.71	0.22	0.015	0.017
Replaces Stryker® 2296-3-212, 2296-003-212					
					Thin
<b>KM3-225</b>	mm	31	9	0.38	0.43
	in	1.22	0.35	0.015	0.017
Replaces Stryker® 2296-3-225, 2296-003-225					
					Thin
<b>KM3-254</b>	mm	41	19.5	0.38	0.68
	in	1.61	0.77	0.015	0.027
Replaces Stryker® 2296-3-254, 2296-003-254					
					Thin
<b>KM3-410</b>	mm	11.5	5.5	0.38	0.43
	in	0.45	0.22	0.015	0.017
Replaces Stryker® 2296-3-410, 2296-003-410					
					Thin

Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
 SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2

**STRYKER®** |

# BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements









**STRYKER®**

Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2

**Compatible with Stryker® Oscillating & Sagittal Saw Systems:**







**SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2**

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM3-411</b>	mm	11.5	9	0.38	0.43
	in	0.45	0.35	0.015	0.017
Replaces Stryker® 2296-3-411, 2296-003-411					
					Thin
<b>KM3-412</b>	mm	18	5.5	0.38	0.43
	in	0.71	0.22	0.015	0.017
Replaces Stryker® 2296-3-412, 2296-003-412					
					Thin
<b>KM3-413</b>	mm	18	2.5	0.38	0.43
	in	0.71	0.10	0.015	0.017
Replaces Stryker® 2296-3-413, 2296-003-413					
					Thin
<b>KM3-414</b>	mm	25	5.5	0.38	0.43
	in	0.98	0.22	0.015	0.017
Replaces Stryker® 2296-3-414, 2296-003-414					
					Thin
<b>KM3-511</b>	mm	25	9	0.38	0.64
	in	0.98	0.35	0.015	0.025
Replaces Stryker® 2296-3-511, 2296-003-511					
					Thin Aggressive Cut
<b>KM3-525</b>	mm	31	9	0.38	0.64
	in	1.22	0.35	0.015	0.025
Replaces Stryker® 2296-3-525, 2296-003-525					
					Thin Aggressive Cut



# BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements

Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
 SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM33-101</b>	mm	10	9	0.51	0.56
	in	0.39	0.35	0.020	0.022
Replaces Stryker® 2296-33-101, 2296-033-101					
					Thick
<b>KM33-102</b>	mm	18.5	16.5	0.51	0.56
	in	0.73	0.65	0.020	0.022
Replaces Stryker® 2296-33-102, 2296-033-102					
					Thick
<b>KM33-103</b>	mm	15	7	0.51	0.56
	in	0.59	0.28	0.020	0.022
Replaces Stryker® 2296-33-103, 2296-033-103					
					Thick
<b>KM33-104</b>	mm	34.5	16.5	0.51	0.56
	in	1.36	0.65	0.020	0.022
Replaces Stryker® 2296-33-104, 2296-033-104					
					Thick
<b>KM33-105</b>	mm	18.5	9	0.51	0.56
	in	0.73	0.35	0.020	0.022
Replaces Stryker® 2296-33-105, 2296-033-105					
					Thick
<b>KM33-106</b>	mm	34.5	13	0.51	0.56
	in	1.36	0.51	0.020	0.022
Replaces Stryker® 2296-33-106, 2296-033-106					
					Thick

Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
 SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2

**STRYKER®**

# BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements









**STRYKER®**

Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2

**Compatible with Stryker® Oscillating & Sagittal Saw Systems:**







**SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2**

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM33-108</b>	mm	39	13	0.51	0.56
	in	1.54	0.51	0.020	0.022
Replaces Stryker® 2296-33-108, 2296-033-108					
					Thick
<b>KM33-111</b>	mm	25	9	0.51	0.56
	in	0.98	0.35	0.020	0.022
Replaces Stryker® 2296-33-111, 2296-033-111					
					Thick
<b>KM33-114</b>	mm	18.5	7	0.51	0.56
	in	0.73	0.28	0.020	0.022
Replaces Stryker® 2296-33-114, 2296-033-114					
					Thick
<b>KM33-115</b>	mm	29.5	7	0.51	0.56
	in	1.16	0.28	0.020	0.022
Replaces Stryker® 2296-33-115, 2296-033-115					
					Thick
<b>KM33-125</b>	mm	31	9	0.51	0.56
	in	1.22	0.35	0.020	0.022
Replaces Stryker® 2296-33-125, 2296-033-125					
					Thick
<b>KM33-206</b>	mm	34.5	12	0.51	0.56
	in	1.36	0.47	0.020	0.022
Replaces Stryker® 2296-33-206, 2296-033-206					
					Thick



# BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements

**Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2**

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM33-212</b>	mm	18	5.5	0.51	0.56
	in	0.71	0.22	0.020	0.022
Replaces Stryker® 2296-33-212, 2296-033-212					
					Thick
<b>KM33-225</b>	mm	31	9	0.51	0.56
	in	1.22	0.35	0.020	0.022
Replaces Stryker® 2296-33-225, 2296-033-225					
					Thick
<b>KM33-410</b>	mm	11.5	5.5	0.51	0.56
	in	0.45	0.22	0.020	0.022
Replaces Stryker® 2296-33-410, 2296-033-410					
					Thick
<b>KM33-411</b>	mm	11.5	9	0.51	0.56
	in	0.45	0.35	0.020	0.022
Replaces Stryker® 2296-33-411, 2296-033-411					
					Thick
<b>KM33-412</b>	mm	18	5.5	0.51	0.56
	in	0.71	0.22	0.020	0.022
Replaces Stryker® 2296-33-412, 2296-033-412					
					Thick
<b>KM33-413</b>	mm	18	2.5	0.51	0.56
	in	0.71	0.10	0.020	0.022
Replaces Stryker® 2296-33-413, 2296-033-413					
					Thick



Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2

**STRYKER®**

# BUSA® Small Bone Sagittal & Oscillating Saw Blades Stryker® Replacements



**Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2**

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM33-414</b>	mm	25	5.5	0.51	0.56
	in	0.98	0.22	0.020	0.022
Replaces Stryker® 2296-33-414, 2296-033-414					
					Thick
<b>KM33-521</b>	mm	25	9	0.51	0.74
	in	0.98	0.35	0.020	0.029
Replaces Stryker® 2296-33-521, 2296-033-521					
					Thick Aggressive Cut

**STRYKER®**

Compatible with Stryker® Oscillating & Sagittal Saw Systems:  
SABO®, Cordless Driver (Sagittal Saw Attachment), RemB™, CORE, TPS™, F1, TPX & Command 2






# BUSA® Small Bone Reciprocating Saw Blades

## Stryker® Replacements

Compatible with Stryker® Reciprocating Saw Systems:  
RemB™, CORE, TPS™, F1, TPX & Command 2

### Stryker®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM5100-137-123</b>	mm	27.0	6.4	0.38	0.51
	in	1.06	0.25	0.015	0.020
Replaces Stryker® 5100-137-123					
					Thin
<b>KM5100-137-133</b>	mm	21.8	6.4	0.38	0.51
	in	0.86	0.25	0.015	0.020
Replaces Stryker® 5100-137-133					
					Thin
<b>KM5100-137-233</b>	mm	27.4	6	0.38	0.51
	in	1.08	0.24	0.015	0.020
Replaces Stryker® 5100-137-233					
					Thin





Compatible with Stryker® Reciprocating Saw Systems:  
RemB™, CORE, TPS™, F1, TPX & Command 2

**STRYKER**®

# BUSA® Small Bone Reciprocating Rasps

## Stryker® Replacements

Compatible with Stryker® Reciprocating Saw Systems:  
RemB™, CORE, TPS™, F1, TPX & Command 2


Stryker®			
Order Number	Unit	Length	Max. Cut Width
<b>KM5100-37-113</b>	mm	11	5.5
	in	0.43	0.22
Replaces Stryker® 5100-37-113, 5100-037-113			
			Small Tear Cross Cut
<b>KM5100-37-114</b>	mm	13	6.5
	in	0.51	0.26
Replaces Stryker® 5100-37-114, 5100-037-114			
			Large Tear Cross Cut
<b>KM5100-37-115</b>	mm	11	5.5
	in	0.43	0.22
Replaces Stryker® 5100-37-115, 5100-037-115			
			Small Tear
<b>KM5100-37-116</b>	mm	13	6.5
	in	0.51	0.26
Replaces Stryker® 5100-37-116, 5100-037-116			
			Large Tear

**STRYKER®**

Compatible with Stryker® Reciprocating Saw Systems:  
RemB™, CORE, TPS™, F1, TPX & Command 2

# BUSA® Small Bone Reciprocating Saw Blades Stryker® Replacements

Compatible with Stryker® Reciprocating Saw System: Micro

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-062</b>	mm	27.2	6	0.64	1.00
	in	1.07	0.24	0.025	0.039
Replaces Stryker® 1675-133					
					Thick




# BUSA® Small Bone Reciprocating Rasps

## Stryker® Replacements



STRYKER® | Compatible with Stryker® Reciprocating Saw System: Micro






Compatible with Stryker® Reciprocating Saw System: Micro

Stryker®			
Order Number	Unit	Length	Max. Cut Width
<b>KM-800</b>	mm	13	6.5
	in	0.51	0.26
Replaces Stryker® 1675-116			
			Large Tear
<b>KM-801</b>	mm	11	5.5
	in	0.43	0.22
Replaces Stryker® 1675-115			
			Small Tear
<b>KM-802</b>	mm	11	5.5
	in	0.43	0.22
Replaces Stryker® 1675-113			
			Small Tear Cross-Cut



# BUSA® Small Bone Sagittal Saw Blades ConMed® Linatec® Replacements

Compatible with ConMed® Linatec® Sagittal Saw Systems:  
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

ConMed® Linatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3100C</b>	mm	41.5	14	0.40	0.61
	in	1.63	0.55	0.016	0.024
Replaces ConMed® Linatec® 5023-33, 5023-033, 5023-133, 5023-233, 5053-33, 5053-033, 5053-133, 5053-233					
					Coarse Straight Edge
<b>KM-3101</b>	mm	25.7	9.5	0.40	0.60
	in	1.01	0.37	0.016	0.024
Replaces ConMed® Linatec® 5023-38, 5023-038, 5023-138, 5023-238, 5053-38, 5053-038, 5053-138, 5053-238, Z338					
					Fine Straight Edge
<b>KM-3101C</b>	mm	25.7	9.5	0.38	0.61
	in	1.01	0.37	0.015	0.024
Replaces ConMed® Linatec® 5023-43, 5023-043, 5023-143, 5023-243, 5053-43, 5053-043, 5053-143, 5053-243					
					Coarse Straight Edge
<b>KM-3102</b>	mm	25.7	4.5	0.40	0.60
	in	1.01	0.18	0.016	0.024
Replaces ConMed® Linatec® 5023-39, 5023-039, 5023-139, 5023-239, 5053-39, 5053-039, 5053-139, 5053-239, Z339					
					Fine Straight Edge
<b>KM-3104</b>	mm	16.1	9.5	0.40	0.63
	in	0.63	0.37	0.016	0.025
Replaces ConMed® Linatec® 5044-145, Z332					
					Fine Straight Edge

Compatible with ConMed® Linatec® Sagittal Saw Systems:  
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

**CONMED® LINATEC®**

# BUSA® Small Bone Sagittal Saw Blades

## ConMed® Linvatec® Replacements








**CONMED® LINVATEC®**

Compatible with ConMed® Linvatec® Sagittal Saw Systems:  
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

**Compatible with ConMed® Linvatec® Sagittal Saw Systems:**

**MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems**

**ConMed® Linvatec®**

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3105</b>	mm	16.1	4.5	0.40	0.63
	in	0.63	0.18	0.016	0.025
Replaces ConMed® Linvatec® 5023-044, 5023-144, 5023-244, 5053-044, 5053-144, 5053-244					
					Fine Straight Edge
<b>KM-3106</b>	mm	25.7	5.5	0.40	0.61
	in	1.01	0.22	0.016	0.024
Replaces ConMed® Linvatec® 5023-32, 5023-032, 5023-132, 5023-232, 5053-32, 5053-032, 5053-132, 5053-232					
					Coarse Straight Edge
<b>KM-3107</b>	mm	34.3	6	0.40	0.60
	in	1.35	0.24	0.016	0.024
Replaces ConMed® Linvatec® 5023-41, 5023-041, 5023-141, 5023-241, 5053-41, 5053-041, 5053-141, 5053-241, Z341					
					Fine Straight Edge
<b>KM-3183</b>	mm	25.7	9.5	0.40	0.60
	in	1.01	0.37	0.016	0.024
Replaces ConMed® Linvatec® 5023-183					
					Fine Straight Edge
<b>KM-3231</b>	mm	41.5	19.5	0.40	0.61
	in	1.63	0.77	0.016	0.024
Replaces ConMed® Linvatec® 5023-31, 5023-031, 5023-131, 5023-231, 5053-31, 5053-031, 5053-131, 5053-231, Z331					
					Coarse Straight Edge




# BUSA® Small Bone Sagittal Saw Blades ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Sagittal Saw Systems:  
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

## ConMed® Linvatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3237</b>	mm	25.7	14	0.40	0.61
	in	1.01	0.55	0.016	0.024


Replaces ConMed® Linvatec® 5023-37, 5023-037, 5023-137, 5023-237, 5053-37, 5053-037, 5053-137, 5053-237



Coarse  
Straight Edge

<b>KM-3240</b>	mm	41.5	9.5	0.40	0.60
	in	1.63	0.37	0.016	0.024

Replaces ConMed® Linvatec® 5023-40, 5023-040, 5023-140, 5023-240, 5053-40, 5053-040, 5053-140, 5053-240, Z340



Fine  
Tapered  
Straight Edge

Compatible with ConMed® Linvatec® Sagittal Saw Systems:  
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
& Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

**CONMED® LINVATEC®**

# BUSA® Small Bone Oscillating Saw Blades

## ConMed® Linvatec® Replacements



**CONMED® LINVATEC®**

Compatible with ConMed® Linvatec® Oscillating Saw Systems:  
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Oscillating Saw Systems

**Compatible with ConMed® Linvatec® Oscillating Saw Systems:  
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Oscillating Saw Systems**

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3000</b>	mm	18.5	5.5	0.40	0.60
	in	0.73	0.22	0.016	0.024
Replaces ConMed® Linvatec® 5023-62, 5023-062, 5023-162, 5023-262, 5053-62, 5053-062, 5053-162, 5053-262, Z362					
					Fine Straight Edge
<b>KM-3001</b>	mm	28.5	5.5	0.40	0.60
	in	1.12	0.22	0.016	0.024
Replaces ConMed® Linvatec® 5023-64, 5023-064, 5023-164, 5023-264, 5053-64, 5053-064, 5053-164, 5053-264, Z364					
					Fine Straight Edge
<b>KM-3002</b>	mm	34.7	10	0.40	0.60
	in	1.37	0.39	0.016	0.024
Replaces ConMed® Linvatec® 5023-69, 5023-069, 5023-169, 5023-269, 5053-73, 5053-073, 5053-173, 5053-273, Z373					
					Fine Straight Edge
<b>KM-3003</b>	mm	22	5.5	0.40	0.60
	in	0.87	0.22	0.016	0.024
Replaces ConMed® Linvatec® N/A					
					Fine Straight Edge
<b>KM-3004</b>	mm	27.9	9.5	0.40	0.60
	in	1.10	0.37	0.016	0.024
Replaces ConMed® Linvatec® Z365					
					Fine Straight Edge
<b>KM-3005</b>	mm	18.5	10	0.40	0.60
	in	0.73	0.39	0.016	0.024
Replaces ConMed® Linvatec® 5023-63, 5023-063, 5023-163, 5023-263, 5053-63, 5053-063, 5053-163, 5053-263					
					Fine Straight Edge









## BUSA® Small Bone Oscillating Saw Blades ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Oscillating Saw Systems:  
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
& Smith & Nephew® Small Bone/Micro Oscillating Saw Systems

### ConMed® Linvatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3006</b>	mm	12	5.5	0.40	0.60
	in	0.47	0.22	0.016	0.024
Replaces ConMed® Linvatec® 5023-60, 5023-060, 5023-160, 5023-260, 5053-60, 5053-060, 5053-160, 5053-260, Z360					
					Fine Straight Edge
<b>KM-3008</b>	mm	28.5	18.5	0.40	0.60
	in	1.12	0.73	0.016	0.024
Replaces ConMed® Linvatec® 5023-66, 5023-066, 5023-166, 5023-266, 5053-66, 5053-066, 5053-166, 5053-266, Z366					
					Fine Radial Edge
<b>KM-3261</b>	mm	12	10	0.40	0.60
	in	0.47	0.39	0.016	0.024
Replaces ConMed® Linvatec® 5023-61, 5023-061, 5023-161, 5023-261, 5053-61, 5053-061, 5053-161, 5053-261, Z361					
					Fine Straight Edge
<b>KM-3265</b>	mm	28.5	10	0.40	0.60
	in	1.12	0.39	0.016	0.024
Replaces ConMed® Linvatec® 5023-65, 5023-065, 5023-165, 5023-265, 5053-65, 5053-065, 5053-165, 5053-265					
					Fine Radial Edge

Compatible with ConMed® Linvatec® Oscillating Saw Systems:  
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
& Smith & Nephew® Small Bone/Micro Oscillating Saw Systems

**CONMED® LINVATEC®**

# BUSA® Small Bone Oscillating Saw Blades ConMed® Livatec® Replacements



**CONMED® LIVATEC®**

Compatible with ConMed® Livatec® Oscillating Saw Systems:  
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
& Smith & Nephew® Small Bone/Micro Oscillating Saw Systems

Compatible with ConMed® Livatec® Oscillating Saw Systems:  
MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
& Smith & Nephew® Small Bone/Micro Oscillating Saw Systems




## ConMed® Livatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3170</b>	mm	12	9	0.40	0.60
	in	0.47	0.36	0.016	0.024
Replaces ConMed® Livatec® 5023-70, 5023-170, 5023-270, 5053-70, 5053-170, 5053-270, Z370					
					Fine Crescentic
<b>KM-3171</b>	mm	25.5	13.8	0.40	0.60
	in	1.00	0.55	0.016	0.024
Replaces ConMed® Livatec® 5023-71, 5023-171, 5023-271, 5053-71, 5053-171, 5053-271, Z371					
					Fine Crescentic
<b>KM-3172</b>	mm	32	18	0.40	0.60
	in	1.26	0.71	0.016	0.024
Replaces ConMed® Livatec® 5023-72, 5023-172, 5023-272, 5053-72, 5053-172, 5053-272, Z372					
					Fine Crescentic
<b>KM-3176</b>	mm	32	13.8	0.40	0.60
	in	1.26	0.55	0.016	0.024
Replaces ConMed® Livatec® 5023-76, 5023-176, 5023-276, 5053-76, 5053-176, 5053-276					
					Fine Crescentic

# BUSA® Small Bone Reciprocating Saw Blades

## ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Reciprocating Saw Systems:  
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Reciprocating Saw Systems

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-060</b>	mm	25.5	7.6	0.40	0.64
	in	1.00	0.30	0.016	0.025
Replaces ConMed® Linvatec® 5053-20, 5053-020, 5053-120, 5053-220, Z320					
					Flat Shank
<b>KM-070</b>	mm	38	7.6	0.40	0.64
	in	1.50	0.30	0.016	0.025
Replaces ConMed® Linvatec® 5053-21, 5053-021, 5053-121, 5053-221, Z321					
					Flat Shank
<b>KM-54</b>	mm	25.0	7.4	0.40	0.60
	in	0.984	0.29	0.016	0.024
Replaces ConMed® Linvatec® 5053-54, 5053-054, 5053-154, 5053-254, Z354					
					Round Shank





Compatible with ConMed® Linvatec® Reciprocating Saw Systems:  
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Reciprocating Saw Systems

**CONMED® LINVATEC®**

# BUSA® Small Bone Reciprocating Rasps

## ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Reciprocating Saw Systems:  
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Reciprocating Saw Systems

ConMed® Linvatec®			
Order Number	Unit	Length	Max. Cut Width
<b>KM-80</b>	mm	11	5.5
	in	0.43	0.22
Replaces ConMed® Linvatec® 5053-080, 5053-80			
			Small Tear Round Shank
<b>KM-81</b>	mm	13	6.9
	in	0.51	0.27
Replaces ConMed® Linvatec® 5053-081, 5053-81			
			Large Tear Round Shank
<b>KM-82</b>	mm	11	5.5
	in	0.43	0.22
Replaces ConMed® Linvatec® 5053-082, 5053-82			
			Small Tear Cross Cut Round Shank
<b>KM-83</b>	mm	13	6.9
	in	0.51	0.27
Replaces ConMed® Linvatec® 5053-083, 5053-83			
			Large Tear Cross Cut Round Shank

**CONMED® LINVATEC®**





Compatible with ConMed® Linvatec® Reciprocating Saw Systems:  
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®, BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Reciprocating Saw Systems



# BUSA® Small Bone Sagittal Saw Blades ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Mini-Driver™  
K100 & K200 Sagittal Saw Systems - K120 Attachment

## ConMed® Linvatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-21K</b>	mm	15.1	6.4	0.71	0.89
	in	0.59	0.25	0.028	0.035
Replaces ConMed® Linvatec® K121					
					Fine
<b>KM-22K</b>	mm	31	6.4	0.71	0.89
	in	1.22	0.25	0.028	0.035
Replaces ConMed® Linvatec® K122					
					Fine
<b>KM-23K</b>	mm	15.1	10.4	0.71	0.89
	in	0.59	0.41	0.028	0.035
Replaces ConMed® Linvatec® K123					
					Fine
<b>KM-24K</b>	mm	31	10.4	0.71	0.89
	in	1.22	0.41	0.028	0.035
Replaces ConMed® Linvatec® K124					
					Fine

Compatible with ConMed® Linvatec® Mini-Driver™  
K100 & K200 Sagittal Saw Systems - K120 Attachment




**CONMED® LINVATEC®**

# BUSA® Small Bone Sagittal Saw Blades ConMed® Linvatec® Replacements



**CONMED® LINVATEC®** | Compatible with ConMed® Linvatec® Mini-Driver™  
 K100 & K200 Sagittal Saw Systems - K120 Attachment






Compatible with ConMed® Linvatec® Mini-Driver™  
 K100 & K200 Sagittal Saw Systems - K120 Attachment

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-31K</b>	mm	15.1	6.4	0.38	0.56
	in	0.59	0.25	0.015	0.022
Replaces ConMed® Linvatec® K131					
					Fine
<b>KM-33K</b>	mm	15.1	10.4	0.38	0.56
	in	0.59	0.41	0.015	0.022
Replaces ConMed® Linvatec® K133					
					Fine
<b>KM-34K</b>	mm	31	10.4	0.38	0.56
	in	1.22	0.41	0.015	0.022
Replaces ConMed® Linvatec® K134					
					Fine



# BUSA® Small Bone Sagittal Saw Blades ConMed® Livatec® Replacements

Compatible with ConMed® Livatec® Mini-Driver™  
K100 & K200 Sagittal Saw Systems - K220 Attachment

ConMed® Livatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-41K</b>	mm	14.2	6.4	0.71	0.97
	in	0.56	0.25	0.028	0.038
Replaces ConMed® Livatec® K141					
					Fine
<b>KM-43K</b>	mm	14.2	10.2	0.71	0.97
	in	0.56	0.40	0.028	0.038
Replaces ConMed® Livatec® K143					
					Fine
<b>KM-44K</b>	mm	30	10.2	0.71	0.97
	in	1.18	0.40	0.028	0.038
Replaces ConMed® Livatec® K144					
					Fine
<b>KM-48K</b>	mm	27.4	9.5	0.38	0.65
	in	1.08	0.37	0.015	0.026
Replaces ConMed® Livatec® K148					
					Fine
<b>KM-51K</b>	mm	14.2	6.4	0.38	0.65
	in	0.56	0.25	0.015	0.026
Replaces ConMed® Livatec® K151					
					Fine







Compatible with ConMed® Livatec® Mini-Driver™  
K100 & K200 Sagittal Saw Systems - K220 Attachment

**CONMED® LIVATEC®**

# BUSA® Sternum Saw Blades


## Stryker® Replacements

Compatible with Stryker® Sternum Saw Systems: 8, 7,6, 5, 4, 2000, EHD & OrthoPower

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>ST6207-97-101</b>	mm	30.5	6.5	0.64	1.02
	in	1.20	0.26	0.025	0.040
Replaces Stryker® 6207-97-101, 6207-097-101					
					1 Per Pack Sterile
<b>ST6207-97-101N</b>	mm	30.5	6.5	0.64	1.02
	in	1.20	0.26	0.025	0.040
Replaces Stryker® 6207-97-101, 6207-097-101					
					1 Per Pack Non-Sterile
<b>ST296-97-102</b>	mm	31.5	6.2	0.64	1.01
	in	1.24	0.24	0.025	0.04
Replaces Stryker® 296-91-102, 296-097-102					
					1 Per Pack Sterile
<b>ST296-97-102N</b>	mm	31.5	6.2	0.64	1.01
	in	1.24	0.24	0.025	0.04
Replaces Stryker® 296-91-102, 296-097-102					
					1 Per Pack Non-Sterile
<b>KM-97</b>	mm	32	6.4	0.79	1.07
	in	1.26	0.25	0.031	0.042
Replaces Stryker® 298-97-100, 298-097-100					
					1 Per Pack Sterile
<b>KM-97N</b>	mm	32	6.4	0.79	1.07
	in	1.26	0.25	0.031	0.042
Replaces Stryker® 2296-3-105, 2296-003-105					
					1 Per Pack Non-Sterile



STRYKER® | Compatible with Stryker® Sternum Saw Systems: 8, 7,6, 5, 4, 2000, EHD & OrthoPower





# BUSA® Sternum Saw Blades Stryker® Replacements

Compatible with Stryker® Sternum Saw Systems: 8, 7, 6, 5, 4, 2000, EHD & OrthoPower

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-98</b>	mm	32	6.4	0.64	1.02
	in	1.26	0.25	0.025	0.040
Replaces Stryker® 298-97-101, 298-097-101					
					1 Per Pack Sterile
<b>KM-98N</b>	mm	32	6.4	0.64	1.02
	in	1.26	0.25	0.025	0.040
Replaces Stryker® 298-97-101, 298-097-101					
					1 Per Pack Non-Sterile


**STRYKER®** | Compatible with Stryker® Sternum Saw Systems: 8, 7, 6, 5, 4, 2000, EHD & OrthoPower

# BUSA® Sternum Saw Blades Stryker® Replacements



**STRYKER®** | Compatible with Stryker® Sagittal Saw Systems: 8, 7, 6, 5, 4, 2000 & EHD





Compatible with Stryker® Sagittal Saw Systems: 8, 7, 6, 5, 4, 2000 & EHD

Stryker®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>BR2108-125</b>	mm	45	42	0.64	1.02
	in	1.77	1.65	0.025	0.040
Replaces Stryker® 2108-125, 2108-125-000					
					
<p>1 Per Pack Sterile Tooth Offset Repeat Sternotomy Blade</p>					

# BUSA® Sternum Saw Blades

## ConMed® Linvatec® Replacements

Compatible with ConMed® Linvatec® Sternum Saw Systems

ConMed® Linvatec®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-31</b>	mm	54	9.9	0.64	0.97
	in	2.13	0.39	0.025	0.038
Replaces ConMed® Linvatec® 5059-31, 5059-031, 5059-531, 5059-031-01					
					Sterile
<b>KM-31N</b>	mm	54	9.9	0.64	0.97
	in	2.13	0.39	0.025	0.038
Replaces ConMed® Linvatec® 5059-31, 5059-031, 5059-531, 5059-031-01					
					Non-Sterile
<b>KM-32</b>	mm	32.8	9.4	0.64	0.99
	in	1.29	0.37	0.025	0.039
Replaces ConMed® Linvatec® 5059-32, 5059-032, 5059-532, 5059-032-01					
					Sterile
<b>KM-32N</b>	mm	32.8	9.4	0.64	0.99
	in	1.29	0.37	0.025	0.039
Replaces ConMed® Linvatec® 5059-32, 5059-032, 5059-532, 5059-032-01					
					Non-Sterile

Compatible with ConMed® Linvatec® Sternum Saw Systems

CONMED® LINVATEC®

# BUSA® Sternum Saw Blades

## ConMed® Linvatec® Replacement & Proprietary



**CONMED® LINVATEC®**


Compatible with ConMed® Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair;  
 ConMed® Linvatec® Micro 100™ Sagittal Saw System; ConMed® Linvatec® Repeat Sternotomy Saw System;  
 BUSA®, MicroAire® & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

**Compatible with ConMed® Linvatec® Oscillating Saw Systems: VersiPower®, Series 4, Series 3 & Orthair; ConMed® Linvatec® Micro 100™ Sagittal Saw System; ConMed® Linvatec® Repeat Sternotomy Saw System; BUSA®, MicroAire® & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems**

**ConMed® Linvatec®**

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-279</b>	mm	39.5	55	0.40	0.79
	in	1.56	2.17	0.016	0.031


Replaces ConMed® Linvatec® 5059-279, Stryker® 5301-10-279, DePuy® Synthes® 5608-13, MicroAire® Z0-079



Sterile  
Repeat Sternotomy Blade  
Opposing Teeth

<b>KM-279PZ</b>	mm	26	39	0.40	0.71
	in	1.02	1.54	0.016	0.028

Replaces ConMed® Linvatec® N/A



Sterile  
Repeat Sternotomy Blade



# BUSA® Sternum Saw Blades


## Proprietary

Compatible with ConMed® Linvatec® Small Bone/Micro Sagittal Saw Systems:  
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

### ConMed® Linvatec®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-3179SN</b>	mm	42.5	54	0.40	0.69
	in	1.67	2.13	0.016	0.027


Replaces ConMed® Linvatec® N/A



Non-Sterile  
 Coarse  
 Opposing Teeth

<b>KM-3231SN</b>	mm	41.9	19	0.40	0.71
	in	1.65	0.75	0.016	0.028


Replaces ConMed® Linvatec® N/A



Non-Sterile

<b>KM-3231SS</b>	mm	41.9	19	0.40	0.71
	in	1.65	0.75	0.016	0.028

Replaces ConMed® Linvatec® N/A




Sterile

Compatible with ConMed® Linvatec® Small Bone/Micro Sagittal Saw Systems:  
 MicroPower®, MicroChoice®, Micro 100™ & Micro E®; BUSA®, MicroAire®  
 & Smith & Nephew® Small Bone/Micro Sagittal Saw Systems

**CONMED® LINVATEC®**

# BUSA® Sternum Saw Blades Codman® (Symmetry) Replacements

Compatible with Codman® Sternum Saw Systems



Codman®					
Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-7010</b>	mm	34	9.4	0.89	1.20
	in	1.34	0.37	0.035	0.047
Replaces Codman® 26-7010					
					1 Per Pack Sterile

**CODMAN®** | Compatible with Codman® Sternum Saw Systems

# BUSA® Sternum Saw Blades Replacements Terumo®

Compatible with Sarns® Sternal Saw II

## Codman®

Order Number	Unit	Length	Width	Material Thickness	Cut Thickness
<b>KM-278</b>	mm	32	6.4	0.79	1.02
	in	1.26	0.25	0.031	0.040
Replaces Terumo® 5589					
					1 Per Pack Sterile
<b>KM-278N</b>	mm	32	6.4	0.79	1.02
	in	1.26	0.25	0.031	0.040
Replaces Terumo® 5756					
					1 Per Pack Non-Sterile